

Sheet No.

**GT200-PE022E** Oil

## Determination of acid number in diesel engine oil 1/6

\*This application sheet is provided as reference, and does not assure the measurement results. Please consider analysis environment, external factors and sample nature for optimal conditions before the measurement.

### Outline

Acid number in diesel engine oil is determined with titration by potassium hydroxide in 2-propanol titrant after dissolving new or used oil into titration solvent contains toluene, 2-propanol and small amount of water. The titration result is used as reference of oxidation and deterioration state of the oils for example.

Titration Type : Non-aqueous Neutralization, Titration mode: INF, Detection: pH/mV  
 ◆Reference : **ASTM D664-07** Standard Test Method for Acid Number of Petroleum Products by Potentiometric Titration

### Apparatus

Automatic titrator : GT-200  
 Electrodes : Reference Electrode sleeve type, Glass electrode  
 Reference electrode solution : 3mol/L, Lithium chloride in ethanol  
 Buret size : 10ml

### Reagents

[ Titration solvent ]

■0.1mol/L- Potassium hydroxide in 2-propanol for testing neutralization number in oil

[ Prepared reagents ]

■Titration solvent : mixed 500ml of toluene, 495ml of 2-propanol and 5ml of pure water

■3mol/L of Lithium chloride in ethanol : Dissolve 12.7g of lithium chloride, special grade reagent, in ethanol, special grade reagent, and dilute the solution to 100ml by the ethanol.

### Analytical Procedure

[ Blank measurement ]

- (1) Add 125ml of the titration solvent into a 200ml beaker by a measuring cylinder.
- (2) Titrate with 0.1mol/L . potassium hydroxide in 2-propanol titrant

[ Sample measurement ]

- (1) Add proper size of sample decided by the method depending on the acid number of the sample into a 200ml beaker. 5g +/-0.5g in this sample.
- (2) Add 125ml of the titration solvent into the above mentioned 200ml beaker by a measuring cylinder.
- (3) Titrate with 0.1mol/L . potassium hydroxide in 2-propanol titrant

Sheet No.

**GT200-PE022E Determination of acid number in diesel engine oil** \_\_\_\_\_ 2/6

[ Calculation ]

**Acid number ( mgKOH/g ) = ( A1 – BL ) × M × E × f × FW / S × R**
**(Used prefixed formula on GT-200)**

A1 : Titration volume of 0.1mol/L- potassium hydroxide in 2-propanol titrant for sample measurement (ml)

BL : Titration volume of 0.1mol/L- potassium hydroxide in 2-propanol titrant for blank measurement (ml)

M : Molarity of 0.1mol/L- potassium hydroxide in 2-propanol titrant (0.1)

E : Equivalent number of 0.1mol/L- potassium hydroxide in 2-propanol titrant (1)

F : Factor of 0.1mol/L- potassium hydroxide in 2-propanol titrant

FW : Formula weight of potassium hydroxide (56.1)

S : Sample size(g)

R : Dilution rate (1)

### Other Requirements

- When measuring samples with pH detection, calibrate the apparatus by three standards, pH 7, 4 and 11 before measurement. Select “Sleeve type liquid: 3.3M KCL (GTRS10B)” and “Three point calibration (Input pH)” on the “pH Calibration” of GT-200.
- For using 10ml Buret, set the volume by “Setting” on the Automatic Buret’ s software.
- After a measurement, wash the electrodes by the titration solvent and immerse them in pure water for 5min. as conditioning.
- Confirm reagent labels and safety data sheets for safety
- Wear protective equipment (eye protector, gloves and others.) when handling reagents.

### Measurement Results

Detection : mV

	Sample size(g)	Titrant (ml)	Results(mg KOH/g)
1	5.0300	2.7613	2.93
2	5.0490	2.7581	2.91
3	5.0464	2.7849	2.95

N                3  
 Average        2.93  
 SD                0.016  
 RSD(%)        0.54

Detection : pH

	Sample size(g)	Titrant (ml)	Results(mg KOH/g)
1	5.0175	2.5524	2.91
2	5.0230	2.5016	2.88
3	5.0049	2.6902	2.89

N                3  
 Average        2.90  
 SD                0.017  
 RSD(%)        0.58

Acid number in diesel engine oil (10W-30) is measured by GT-200.

Average of three measurements is around 2.9mgKOH / g. The results are repeatable on both mV and pH detections.

Sheet No.


**GT200-PE022E** Determination of acid number in diesel engine oil 3/6

ID No. : 7 GT No.1

User : GT-200

Measurement : 2014/08/21 15:00  
 Sample Name : BLANK

Type : Sample Titr  
 Sample Size : 125 [ml]



C1 : 0.1357 [ml]

A1 : 0.1357 [ml] -219 [mV]

P-initial : -17 [mV]  
 Start : 0 [ml] -17 [mV]  
 End : 0.45 [ml] -470 [mV] Measuring Time : 4'47"

File No. : 14 OIL / Acid Number  
 Titr File No. : 39 Acid Number / Blank  
 Mode : INF End1, End1 Width : -500 [mV] ± 1000 [mV]  
 Detect : mV1  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 10 [mV]  
 dT : 10 [sec]  
 DL : 600 [mV/ml]  
 DetCnt : 2  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml] C1 : A1  
 [ml]  
 Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1

Buret Injection Speed : 500 [ul/sec]

Sheet No.


**GT200-PE022E** Determination of acid number in diesel engine oil 4/6

ID No. : 8 GT No.1

User : GT-200

Measurement : 2014/08/21 15:28  
 Sample Name : Engine Oil

Type : Sample Titr  
 Sample Size : 5.0300 [g]



C1 : 2.93 [mgKOH/g]

A1 : 2.7613 [ml] -286 [mV]

P-initial : 19 [mV]  
 Start : 0 [ml] 19 [mV]  
 End : 3.864 [ml] -460 [mV] Measuring Time : 13'58"

File No. : 14 OIL / Acid Number  
 Titr File No. : 6 OIL / Acid Number  
 Mode : INF End1, End1 Width : -300 [mV] ± 500 [mV]  
 Detect : mV1  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 10 [mV]  
 dT : 10 [sec]  
 DL : 200 [mV/ml]  
 DetCnt : 20 C1 : (A1-BL)\*M\*E\*f\*FW/S\*R  
 Vmax : 20 [ml] [mgKOH/g]  
 Vover : 0.2 [ml]

Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1 BL : 0.1357 [ml]  
 FW : 56.1 R : 1

Buret Injection Speed : 500 [ul/sec]

Sheet No.


**GT200-PE022E** Determination of acid number in diesel engine oil 5/6

ID No. : 1GT No.1

User : GT-200

Measurement : 2014/08/20 10:10  
 SampleName : Blank

Type : Sample Titr  
 SampleSize : 125 [ml]



C1 : 0.1261 [ml]

A1 : 0.1261 [ml] 11.229 [pH]

P-initial : 6.939 [pH]  
 Start : 0 [ml] 6.939 [pH]  
 End : 1.75 [ml] 14 [pH] Measuring Time: 6'15"

FileNo. : 14 OIL/ Acid Number  
 Titr FileNo. : 40 AcidNumber / BlankpH  
 Mode : INF End1, End1 Width: 11[pH] ± 2 [pH]  
 Detect : pH  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 0.2 [pH]  
 dT : 10 [sec]  
 DL : 20 [pH/ml]  
 DetCnt : 3  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml] C1: A1  
 [ml]

Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1

BuretInjectionSpeed: 500 [ul/sec]

Sheet No.


**GT200-PE022E** Determination of acid number in diesel engine oil 6/6

ID No. : 1 GT No.1

User : GT-200

Measurement : 2014/08/21 10:32  
 Sample Name : Diesel oil

Type : Sample Titr  
 Sample Size : 5.0175 [g]



C1 : 2.91 [mgKOH/g]

A1 : 2.7324 [ml] 11.646 [pH]

P-initial : 5.74 [pH]  
 Start : 0 [ml] 5.74 [pH]  
 End : 3.506 [ml] 14 [pH] Measuring Time : 12'13"

File No. : 14 OIL / Acid Number  
 Titr File No. : 36 OIL / Acid Number pH  
 Mode : INF End1, End1 Width : 11 [pH] ± 2 [pH]  
 Detect : pH  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 0.2 [pH]  
 dT : 10 [sec]  
 DL : 5 [pH/ml]  
 DetCnt : 10  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml]

C1 : (A1-BL)\*M\*E\*f\*FW/S\*R [mgKOH/g]

Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1 BL : 0.1261 [ml]  
 FW : 56.1 R : 1

Buret Injection Speed : 500 [ul/sec]